


OCA 15EC

The entry level contact angle measuring and contour analysis system






analysis of the roll off angle with an OCA 15EC on electronic tilting base unit TBU 95

The **OCA 15EC** is the entry level measuring device for professional contact angle measurements and drop shape analysis. The OCA 15EC can be easily disassembled for transportation in the optional case. For a convenient access to the surface and interfacial measurement techniques two pre-configured packages are available which include all necessary components to start right away.



Contact angle measurement with an OCA 15EC with double direct dosing system DD-DM




liquid temperature control unit TFC 100Pro for temperatures of -10...100 °C

The **OCA 15EC Package 1** consists in addition to the base unit of a single direct dosing system SD-DM, one electronic syringe unit ESr-N and the software module SCA 20.


The **OCA 15EC Package 2** consists in addition to the base unit of a double direct dosing system DD-DM, two electronic syringe units ESr-N and the software module SCA 20.

Main features:


- Sample table horizontally movable (magnetic slide system) and adjustable vertically with precision mechanics
- high performance 6-times zoom lens
- integrated continuous fine focus, and adjustable observation angle
- video measuring system with USB camera (max. 311 frames/s), high-speed camera easily upgradable (up to 3000 frames/s)
- LED-lighting with manual and software controlled intensity including automatic temperature drift compensation




interfacial tension determination with pendant drop method



schematic contour analysis of a pendant drop



contact angle of a sessile drop



SCA 22 — analysis of the interfacial tension

Software for efficient work

The SCA software, designed for Microsoft Windows®, is the modular program for all OCA instruments. The available software modules for the OCA 15EC are:

SCA 20 — contact angle

- video based measurement and presentation of the static and dynamic contact angle on plane, convex, and concave surfaces
- automatic measurement of the contact angle hysteresis
- record/store of image sequences
- statistics and measurement error analysis
- liquids and solids database

SCA 21 — surface free energy

- analysis of the surface free energy of solids as well as its components (e.g. dispersive, polar and hydrogen bond parts, acid and base portions) according to nine different theories
- calculation and representation of wetting envelopes and work of adhesion/contact angle diagrams

SCA 22 — surface/interfacial tension


- analysis of the surface and interfacial tension, as well as their polar and dispersive parts, based on the analysis of the shape of pendant drops

SCA 23 — liquid bridge analysis

- analysis of the surface and interfacial tension based on the evaluation of the lamella contour
- innovative liquid bridge analysis of 3 phase systems

SCA 26 — oscillation / relaxation

- analysis of the real and imaginary part of the interfacial dilatational modulus based on the oscillating or relaxing contour of pendant drops



SCA 20 and 21 — contact angle measurement and surface free energy determination

Technical data

Max. sample dimensions (L x W x H):	• 220 x ∞ x 70 mm ³
Sample table dimensions (L x W):	• 100 x 100 mm ²
Traversing range of sample table in X-Y-Z direction:	• 110 x 90 x 42 mm ³
Measuring range for contact angles:	• 0...180°; ± 0.1° measuring precision of the video system
Measuring range for surface and interfacial tensions:	• 1·10 ⁻² ... 2·10 ³ mN/m; resolution: min. ± 0.01 mN/m
Max. sample weight:	• 3.0 kg; 15.0 kg with clamped sample table
Optics and image processing system:	<ul style="list-style-type: none"> • LED-lighting with manual and software controlled intensity including automatic temperature drift compensation • USB 2.0 camera, max. resolution 752 x 480 pixel, max. frame rate 311 frames/s • 6-fold zoom lens with integrated fine focus (± 6 mm) • field of view: 1.05 x 0.66...6.72 x 4.25 mm² • optical distortion: < 0.05 %
Dimensions (L x B x H):	• 550 x 160 x 365 mm ³
Weight:	• 14 kg
Power supply:	• 100...240 V AC; 50...60 Hz; 70 W

Accessories (excerpt)

manual direct dosing systems **SD-DM** and **DD-DM** with corresponding electronic syringe units **ESr** • electronic tilting base unit **TBU 95** • temperature and environmental control chambers (-30...700 °C) • syringe heating device **SHD** (up to 90 °C) • holders for foils or papers **FSH 30** and **FSC 80/150** • sample table with holding clamps **STC 100** • film or foil sample stage **FHM 100** • suction plate **SP 100** for holding thin flexible samples flat on the stage • electro wetting platform **EWP 100**

For more information about a tailor made solution to your surface chemistry requirements, please contact us. We will be pleased to provide a quotation, obligation free, for your instrument system.

DataPhysics Instruments GmbH • Raiffeisenstraße 34 • 70794 Filderstadt, Germany
 phone +49 (0)711 770556-0 • fax +49 (0)711 770556-99
sales@dataphysics.de • www.dataphysics.de

Representante en España y Portugal

NEURTEK
 instruments

Pol. Ind. Azitain 3A
 20600 Eibar

Tel. 943820082
 Fax 943820157
comercial@neurtek.es
sat@neurtek.es
www.neurtek.com

Delegaciones: Barcelona, Madrid, Sevilla, Vigo y Oporto