


NEURTEK

i n s t r u m e n t s

**HOLDER FOR DOSING NOZZLES**

The holder for dosing nozzles for manual operation can be freely positioned on the machine. Dosage is controlled via the operator panel of the TOPAS M.

**TOPAS M**

The automatic dosing system TOPAS M is universally suitable for use with all manual grinders. Four 500 ml glass bottles with screw caps hold 3 x suspension and 1 x lubricant or 4 x Dia-Complete and can be filled via a separate lid. They are easily removed for cleaning.

The holder for dosing nozzles for manual dosage may be freely positioned on the machine. The dosage interval can be variably set via the control panel.

The backflushing function keeps the hoses clean.

** Operable only with electric equipment*

SUSPENSION CONTAINERS

The four glass bottles with screw cap have a capacity of 500 ml. They can be filled during operation via a separate inlet in the lid.

They are automatically arrested and are easily removed.



BASIC MODULE



TOPAS M



- » Automatic dosing unit for
3x diamond suspension and 1x lubricant or
4x Dia-Complete
- » Dosing interval adjustable
- » Reverse rinsing function
- » Manual or automatic mode
- » For all grinders and polishers
- » Suspension containers refillable and removable

100-240 V / 50/60 Hz (1Ph/N/PE)

Order No.: M6300050

HOLDER FOR DOSING NOZZLES (MANUAL)



Order No.: A6300051

- » 4 nozzles
- » Combinable with SAPHIR grinders and other brands
- » Arm W x H x D: 250 x 55 x 50 mm

SPECIFICATIONS

Bottles	4
Capacity	500 ml
Power supply	0.2 kVA
Connected load	100-240 V / 50/60 Hz (1Ph/N/PE)

WxHxD ~	158 x 686 x 223 mm
Weight	9 kg

